

Contents

| | | |
|----------|---|----------|
| 1 | Necessary Hardware | 2 |
| 2 | Setup | 3 |
| 2.1 | Hardware Setup | 3 |
| 2.1.1 | Synchronize the two USRP1s | 3 |
| 2.1.2 | Attach LFTX and LFRX daughterboards in appropriate mother- board slots | 3 |
| 2.2 | Build "amplifier" circuit for TXE pathway | 4 |
| 2.3 | Connect SMA bulkhead cables to enclosure output/input slots and connect internal feedback cables | 4 |
| 2.4 | Connect USRP1 outputs and inputs to external hardware | 5 |
| 3 | Software Setup | 5 |
| 3.1 | GNU Radio Setup | 5 |
| 3.2 | Downloading and Installing gr-MRI | 5 |
| 3.3 | Editing GNU Radio settings for stock sequences | 6 |
| 4 | Using gr-MRI | 7 |
| 4.1 | Overview | 7 |
| 4.2 | RF Pulse Calibration | 7 |
| 4.3 | Gradient Calibration | 7 |
| 4.4 | Stock Imaging Sequences | 7 |
| 4.4.1 | Slice Selective Gradient Echo: <i>gradecho.py</i> | 7 |
| 4.4.2 | Slice Selective Spin Echo: <i>slicespin.py</i> | 7 |
| 4.4.3 | Slice Selective Inversion Recovery: <i>invrecov.py</i> | 7 |

1 Necessary Hardware

In order to construct the USRP1 based MRI spectrometer as we have, you will need the following hardware:

1. Ettus Research:

- (2x) **Ettus Research USRP1** (\$707 ea.): <http://www.ettus.com/product/details/USRPPKG>

Note: This package comes with a USB cable, power supply, and 2 SMA bulk-head cables

- (4x) **Ettus Research LFRX** (\$76 ea.): <http://www.ettus.com/product/details/LFRX>
- (4x) **Ettus Research LFTX** (\$76 ea.): <http://www.ettus.com/product/details/LFTX>
- (?6x) **SMA Bulkhead Cable** (\$20 ea): <http://www.ettus.com/product/details/SMA-Bulkhead>

2. Digikey (amp circuit):

- **item 1**
- **item 2**
- **item 3**

3. Miscellaneous:

- SMA-SMA cable to connect clocks
- (2x) short SMA-SMA cable for internal loopbacks
- misc electrical wires
- 8x BNC to SMA connectors (female to male)

2 Setup

2.1 Hardware Setup

2.1.1 Synchronize the two USRP1s

The clocks on the different USRP1s may be slightly different, and without synchronization pulses sent out of both boards will shift relative to one another in time. It is necessary to synchronize the on board clocks to enable sequenced operation. The following instructions have been slightly modified from the [GNU Radio clock synchronization instructions](#).

1. Decide which board will be the master, and which will be the slave. The gr-MRI package assumes that the RF board will serve as the master, so you should do the same.
2. For the master clock board:
 - (a) Solder an SMA connector into J2002. This is the master clock output. Be careful when soldering the SMA connector so you don't break the delicate trace from J2002 to R2028.
3. For the slave boards:
 - (a) Solder an SMA connector into J2001. This is the clock input. Be careful when soldering the SMA connector so you don't break the delicate trace from J2001 to C927.
 - (b) Move R2029 to R2030. This disables the onboard clock. R2029/R2030 is a 0-ohm resistor.
 - (c) Move C925 to C926.
 - (d) Remove C924.
 - (e) If you want to chain another USRP device off of this one, you can use J2002 to provide a clock out
4. To synchronize the clocks, simply connect the connectors that were soldered on.

2.1.2 Attach LFTX and LFRX daughterboards in appropriate motherboard slots

Insert the LFTX and LFRX daughterboards into the appropriate slots on the motherboard according to Figure 1 below. The slots are also labeled *TX_A*, *TX_B*, *RX_A*, *RX_B* on the motherboard. Screw the daughterboards so they are secure. **For this guide I will refer**

to ports as *m-board d-board_slot_channel*. So for example, The RX_B channel for the RF board's RX_A slot would be *RF RX_A_B*.




Figure 1: Image of Ettus Research USRP1 motherboard, LFTX, and LFRX daughterboard placement

2.2 Build "amplifier" circuit for TXE pathway

Come back to this guy later

2.3 Connect SMA bulkhead cables to enclosure output/input slots and connect internal feedback cables

1. Connect the following ports to the output/input holes in the USRP1 enclosure. It doesn't matter which holes to connect to, but you may want to write which slot corresponds to which daughterboard or desired signal.
 - RF TX_A_A (*RF Channel 1 Signal Out*)
 - RF RX_A_A (*RF Signal RX*)

- RF Master Clock - connect to DC Slave Clock
 - RF TX_B_A (*Synchronizing Pulse*)
 - DC TX_A_A (*Z Gradient*)
 - DC TX_A_B (*Y Gradient*)
 - DC TX_B_B (*X Gradient*)
 - DC Slave Clock - connect to RF Master Clock
 - DC RX_A_A (*RF sync pulse input*)
2. Connect the following internal feedback loops:
- Connect RF TX_A_B to RF RX_B_A (*Readout Window*)
3. Special Connections:
- Connect DC TX_B_A through TXE amplifier output 1 to output slot (*TXE pulse*)
 - Connect DC TX_B_A through TXE amplifier output 2 to DC RX_B_A (*sync feedback pulse*)

2.4 Connect USRP1 outputs and inputs to external hardware

RF amplifiers, TX/RX coil, Preamplifier, Gradient Amplifiers, etc.

3 Software Setup

3.1 GNU Radio Setup

Download and install GNU Radio and all of its dependencies following the instructions at <https://gnuradio.org/redmine/projects/gnuradio/wiki/UbuntuInstall>. This will take some time.

3.2 Downloading and Installing gr-MRI

To download the current version of the gr-MRI package on Linux machine:

1. Open a new terminal and navigate to the directory in which you want to save gr-MRI.

2. type `git init`
3. type `git remote add origin ssh git@bitbucket.org:wgrissom/gnuradio-mri.git`
4. type `git clone git@bitbucket.org:wgrissom/gnuradio-mri.git`

To add the gr-MRI custom blocks for use in GNU Radio:

1. Open a terminal and navigate to `gnuradio-mri/gr_3.7/gr-MRI`, and type:
2. `mkdir build`
3. `cd build`
4. `cmake ../`
5. `make`
6. `sudo make install`
7. `sudo ldconfig`

3.3 Editing GNU Radio settings for stock sequences

You will need to change some of the settings within each of the GNU Radio flowgraphs to allow them to work with your USRP1s. The files to edit are called *spin_echo.grc*, *gradient_echo.grc*, *calibration.grc*.

- connect the two USRP1 radios that you plan to use for open a terminal window and type `uhd_find_devices` to make sure that your computer can see the devices. The function will return the serial numbers of the devices. Look on the motherboard of each radio to see which serial number corresponds to which radio (RF or DC).
- Open the `.grc` file and find the blocks titled "USRP sink." The topmost block is always the RF radio sink, and the block below it is the DC radio sink for the files in `gr-MRI`.
- Double click the top USRP sink to open the settings, shown in Figure 2, and replace the serial number written in the "Device Address" to the serial number of your RF radio.
- Repeat the same procedure for the bottom USRP Sink




Figure 2: Image of Ettus Research USRP1 motherboard, LFTX, and LFRX daughterboard placement

4 Using gr-MRI

4.1 Overview

4.2 RF Pulse Calibration

4.3 Gradient Calibration

4.4 Stock Imaging Sequences

4.4.1 Slice Selective Gradient Echo: *gradecho.py*

4.4.2 Slice Selective Spin Echo: *slicespin.py*

4.4.3 Slice Selective Inversion Recovery: *invrecov.py*